

[All Products](#) / [Lab and Production](#) / [Optical Cleaning](#) / [Gloves, Brushes, and Swabs](#)

## Cotton-Tipped Swab Applicators



- Lint-Free
- 3" or 6" Wooden Stems
- Bags, Boxes, and Cases Available

## Technical Information

## Products

Title	Stock Number	Price	Buy
3" Length Bag, Cotton-Tipped Swabs	<a href="#">#56-926</a>	€5,75	Contact Us
3" Length Box (10 Bags), Cotton-Tipped Swabs	<a href="#">#56-927</a>	€27,00	In Stock
3" Length Case (10 Boxes), Cotton-Tipped Swabs	<a href="#">#56-928</a>	€150,00	3-4 days
6" Length Bag, Cotton-Tipped Swabs	<a href="#">#56-924</a>	€5,75	3-4 days
6" Length Box (10 Bags), Cotton-Tipped Swabs	<a href="#">#39-944</a>	€27,00	3-4 days
6" Length Case (10 Boxes), Cotton-Tipped Swabs	<a href="#">#56-925</a>	€150,00	Contact Us



Copyright 2023 | Edmund Optics, Ltd Unit 1, Opus Avenue, Nether Poppleton, York, YO26 6BL, UK

Phone: 1-800-363-1992 :

[www.edmundoptics.com](http://www.edmundoptics.com)