

TECHSPEC® 100mm Length, T-Mount Extension Tube



Stock **#65-969** 3-4 DAYS

- 1 + **€84,00**

ADD TO CART

Qty 1-4
€84,00

Qty 5+
€74,50

Volume Pricing
Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Extension Tube

Extension Length (mm):
100.0

Total Length (mm):
103.45

Outer Diameter (mm):

48.0

Note:
Inner threads are female T-Mount (M42 x0.75),
allowing optical components to be mounted with [#58-741](#)

Thread Type:
T-Mount

Clear Aperture CA (mm):
37.0

Function:
Extension Tube (Fixed)

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)

Reach 235:
[Compliant](#)