

[See all 18 Products in Family](#)

TECHSPEC[®]

30mm, 60mm Radius, Metric Goniometer



TECHSPEC[®] Goniometers



+1

Stock #66-524 6 In Stock

-

1

+

€468^{.00}

ADD TO CART

Volume Pricing	
Qty 1-2	€468,00 each
Qty 3+	€408,10 each
Need More?	Request Quote



ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General	
Type:	Metric
Resolution (arcmin):	5
Physical & Mechanical Properties	

Technical drawing of a 3D-printed part with the following dimensions and specifications:

- Overall width: 30
- Overall height: 20
- Top flange width: 24
- Top flange height: 4
- Central hole diameter: $\varnothing 10$
- Material: $MA (1/4-20)$, 2mm deep
- Threaded hole: $8 \times M2 (4-7.5)$, 4mm deep
- Overall length: 20.5



B (mm)	English Stage	Matrix Stage	A (mm)	B (mm)	C (mm)	D (mm)
40	#66-S20	#66-S23	45	23	7.4	7
60	#66-S21	#66-S24	67	45	8.4	9
82	#66-S22	#66-S25	89	67	8.4	9