

In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union

SUBMIT

TECHSPEC® M25.5 X 0.5 to C-Mount Rotating Male Adapter



#66-381: M25.5 x 0.5 to C-Mount Rotating Male Adapter



Stock #66-381 3-4 DAYS

1 €55,00

ADD TO CART

Qty 1+

€55,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Lens Accessory

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:

