

TECHSPEC® Aperture Kit for Achromat Prototyping Kit



Aperture Kit for Achromat Prototyping Kit

Stock **#54-720** 3-4 DAYS



1



€146,10

ADD TO CART

Qty 1-5

€146,10

Qty 6+

€135,90

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Prototyping Kit

Function:
Kit

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
