


[See all 21 Products in Family](#)

## Barrel Storage Cap, 10cc Barrels (100 Count)



Stock #25-100 **20+ In Stock**

⊖ 1 ⊕ €83<sup>00</sup>

**ADD TO CART**

### Volume Pricing

Qty 1+	€83,50 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

### Product Downloads

### Regulatory Compliance

Certificate of Conformance:

[View](#)

## Product Details

- Empty Barrels are Available in Clear, Amber, and Opaque Versions
- Tips are Available in Plastic and Steel
- [Dispensing Guns](#) and [Pre-Loaded Norland Adhesive Barrels](#) are Also Available

Dispenser barrels, tips, and other accessories are available for use with dispensing guns and hand plungers. Dispenser Barrels are available in a clear version for normal fluids, an amber version for UV-curing fluids, and an opaque version for light sensitive fluids. Dispenser tips are available in plastic and steel to ensure varying levels of controllable deposits without leakage.

**Note:** Barrels, tips, and adhesive are sold separately from [dispensers](#). [Pre-loaded Norland adhesive barrels](#) are also available.

---