


[See all 21 Products in Family](#)

Clear, 3cc Barrel (100 count), Dispenser Barrel



Stock #25-095 **2 In Stock**

- 1 + €167⁰⁰

ADD TO CART

Volume Pricing

Qty 1+	€167,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

Regulatory Compliance

Certificate of Conformance:

[View](#)

Product Details

- Empty Barrels are Available in Clear, Amber, and Opaque Versions
- Tips are Available in Plastic and Steel
- [Dispensing Guns](#) and [Pre-Loaded Norland Adhesive Barrels](#) are Also Available

Dispenser barrels, tips, and other accessories are available for use with dispensing guns and hand plungers. Dispenser Barrels are available in a clear version for normal fluids, an amber version for UV-curing fluids, and an opaque version for light sensitive fluids. Dispenser tips are available in plastic and steel to ensure varying levels of controllable deposits without leakage.

Note: Barrels, tips, and adhesive are sold separately from [dispensers](#). [Pre-loaded Norland adhesive barrels](#) are also available.

;