

[See all 70 Products in Family](#)

TECHSPEC[®] DAPI (High Contrast) Filter Cube Set - Nikon



Olympus Filter Cube #86-833



Stock #86-371-NKN

CONTACT US

-

1

+

€1.350⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-9	€1.350,00 each
Qty 10+	€1.215,00 each
Need More?	Request Quote

 Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Type:

Fluorescence Filter Kit

Compatible Fluorophore:

DAPI (High Contrast)

Excitation Filter:

[#84-094](#)

#84-095	Emission Filter:
#86-330	Dichroic Filter:
Nikon	Compatible Microscope Brand:

Optical Properties

Fused Silica (Corning 7980)	Substrate: □
Hard Coated	Coating:
382 - 393	Excitation Wavelength (nm):
417 - 477	Emission Wavelength (nm):
409.00	Dichroic Cut-On Wavelength (nm):
382 - 477	Wavelength Range (nm):

Regulatory Compliance

View	Certificate of Conformance:
------	-----------------------------

Product Details

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).