

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

Deep Blue 47 Kodak Wratten Color Filter



Stock #53-700 ☒ In Stock



1



€115,00

ADD TO CART

Qty 1-10

€115,00

Qty 11+

€101,00

Volume Pricing

[Request Quote](#)

 Prices shown are exclusive of VAT/local taxes

Get Product Downloads



Dimensions (mm):
75.0 x 75.0

Glass/Filter Number:
47.00

Substrate: ☐
Wratten 2

Length (mm): 75.00
Coating: Uncoated
Color: Deep Blue
Thickness (mm): 0.10
Type: Color Filter
Width (mm): 75.00

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 235: Compliant
