

[See all 1 Products in Family](#)

High Performance Vibration Isolation System

See More by [TMC](#)



Stock **#87-226** [CONTACT US](#)

⊖ 1 ⊕ €3.975⁰⁰

ADD TO CART

Volume Pricing

Qty 1-4	€3.975,00 each
Qty 5+	€3.775,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Universal **Type:**

Physical & Mechanical Properties

3.0 **Thickness (inches):**

20.0 **Length (inches):**

18.0 **Width (inches):**

Construction:
Electropolished stainless steel sides and base
Textured powder coat carbon steel top plate

80 (nitrogen or air) **Feed Pressure (psi):**

3.0 **Height (inches):**

67 (150 lb) **Load Capacity (kg):**

22 (50 lb) **Weight (kg):**

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)

Product Details

- Passive Vibration Isolation
- Lightweight, Compact Profile for Easy Portability
- Self-Leveling Design with Height Control Valves

The TMC™ High Performance Vibration Isolation System is designed to remove the unwanted effects of vibration from challenging microscopy applications. The lightweight, low-profile design weighs only 22kg, easing portability and lab integration. The TMC™ High Performance Vibration Isolation System consists of an electropolished stainless steel chassis and powder finish steel top plate. The self-leveling TMC™ High Performance Vibration Isolation System also includes height control valves for modifying the height of the platform.