

High Performance Vibration Isolation System

See More by [TMC](#)



Stock **#87-226** [CONTACT US](#)

-

1

+

€3.975⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-4	€3.975,00 each
Qty 5+	€3.775,00 each
Need More?	Request Quote

i Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Type:

Universal

Physical & Mechanical Properties

Thickness (inches):

3.0

20.0	Length (inches):
18.0	Width (inches):
Construction: Electropolished stainless steel sides and base Textured powder coat carbon steel top plate	
80 (nitrogen or air)	Feed Pressure (psi):
3.0	Height (inches):
67 (150 lb)	Load Capacity (kg):
22 (50 lb)	Weight (kg):

Regulatory Compliance

Compliant	RoHS 2015:
View	Certificate of Conformance:

Product Details

- Passive Vibration Isolation
- Lightweight, Compact Profile for Easy Portability
- Self-Leveling Design with Height Control Valves

The TMC™ High Performance Vibration Isolation System is designed to remove the unwanted effects of vibration from challenging microscopy applications. The lightweight, low-profile design weighs only 22kg, easing portability and lab integration. The TMC™ High Performance Vibration Isolation System consists of an electropolished stainless steel chassis and powder finish steel top plate. The self-leveling TMC™ High Performance Vibration Isolation System also includes height control valves for modifying the height of the platform.