

InfiniProbe S-25 Video Microscope (0-2.5X)

See More by [Infinity Photo-Optical Company](#)



Stock **#58-806** [CONTACT US](#)

- 1 + **€2.030,00**

ADD TO CART

Qty 1+
€2.030,00

Volume Pricing
Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Horizontal Field of View, 1/2" Sensor:
2.56mm

Horizontal Field of View, 2/3" Sensor:
3.52mm

Weight (g):
158.76

Type:
Long Distance Microscope

Primary Magnification PMAG:

0X - 2.5X
Maximum Sensor Format: 2/3"
Working Distance (mm): 32 - ∞
Mount: C-Mount

Regulatory Compliance

RoHS 2015: Compliant
Reach 224: Compliant
Certificate of Conformance: View
