

InfiniTube™ Ultima, Standard Version

See More by [Infinity Photo-Optical Company](#)



Stock **#34-482** **1 In Stock**

-

1

+

€2.190⁰⁰

ADD TO CART

Volume Pricing	
Qty 1+	€2.190,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Note: This item requires accessories for use | [Learn More](#)

Product Downloads

General	
Tube Lens	Type:
Infinity	Manufacturer:
Objective Design: Nikon CFI60, Olympus UIS/UIS2, Mtuotoyo, Edmund, and Achrovid	

Sensor	
35mm (43.3mm Image Circle)	Sensor Format:
Threading & Mounting	
T-Mount	Mount:
Regulatory Compliance	
Compliant	RoHS 2015:
Compliant	Reach 224:
View	Certificate of Conformance:

Product Details

- Color-Matched for all Infinity Corrected Objectives
- Compatible with Full-Frame 35mm Sensors
- Maintains Perfect Parfocal Performance
- Requires One Objective and Adapter for Use

InfiniTube™ Ultima Microscope Tube Lenses are simplified versions of the Infinity InFocus™ Modules. Except for the InFocus' expanded dynamic focusing, all three color-corrected Ultima modules share the same optics and accessories - flexibly accommodating even the largest formats for increased resolution.

Note: Adapters and accessories not included. One objective and adapter required for use.