

[See all 3 Products in Family](#)

**TECHSPEC® M26/C-Mount Turret**



C-Mount Manual Objective Turrets

Stock #19-534 **9 In Stock**

⊖ 1 ⊕ €900<sup>00</sup>

**ADD TO CART**

Volume Pricing	
Qty 1-2	€900,00 each
Qty 3+	€855,30 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

**Threading & Mounting**

Thread Type:

C-Mount  
Objective: (4) M26 x 0.706

**Regulatory Compliance**

RoHS 2015:


[Compliant](#)

## Product Details

- C-Mount Tube System Compatibility for Creating Modular Microscope System
- Four Objective Positions
- [Mitutoyo](#) and [Nikon](#) Objective Compatibility

C-Mount Manual Objective Turrets are ideal for creating modular microscope systems using standard [C-Mount Tube System Components](#). Edmund Optics' broad range of C-mount Tube System Components offers maximum versatility when designing a microscope system, with direct C-Mount camera compatibility for imaging integration. With additional accessories, these objective turrets can be adapted to [S-Mount](#), [T-Mount](#), [ID1 \(SM1\)](#), and [Multi-Element Tube Systems](#), expanding their potential uses. C-Mount Manual Objective Turrets are offered with M26 x 0.706 threads for [Mitutoyo Objectives](#) and M25 x 0.75 threads for [Nikon Objectives](#). Each turret has four objective positions to maximize imaging system flexibility.

## Technical Information



#19-534