

Motorized XY Microscope Stage Universal Insert

See More by [Zaber™](#)



Stock **#22-632** **2 In Stock**

-

1

+

€765⁰⁰

ADD TO CART

Volume Pricing	
Qty 1+	€765,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General	
Universal	Type:
Zaber Technologies Inc.	Manufacturer:
Note: Compatible with 50mm x 75mm and 25mm x 75mm Slides; Compatible with 35mm and 60mm Diameter Petri Dishes	
Threading & Mounting	

Compliant

RoHS 2015:

View

Reach 247:

Compliant

- 100 x 120mm or 205 x 205mm Travel
- 100N Centered Load Capacity, 12µm Full Travel Accuracy, 2µm Repeatability, and Up to 85mm/s Speed
- Supported by [µManager](#) Open-Source Microscopy Software
- Available with Integrated, 200 Counts per Revolution (CPR) Motor Mounted Encoder

Note: A24-48 VDC universal power supply, data cables for daisy chaining, and computer interface cables (USB or RS-232) are sold separately as accessories. These stages utilize the same accessories as the Zaber™ High Precision Motorized Stage System.

Technical drawing of the front and side views of the 'C' frame. The front view shows a rectangular frame with a central opening. Dimensions include an overall width of 302, an inner width of 125, and an overall height of 236.3. The inner opening has a height of 106.5 and a width of 50. The side view shows a depth of 215.5. The drawing also indicates the use of M6 (x6) and M2 (x4) screws.