

# Olympus Near-Infrared Single Port Tube Lens, 180mm Focal Length

See More by [Olympus](#)



Olympus Compact Tube Lens, 180mm Focal Length

Stock **#75-251** NEW **1 In Stock**

-

1

+

€1.240<sup>00</sup>

ADD TO CART

| Volume Pricing |                               |
|----------------|-------------------------------|
| Qty 1+         | €1.240,00 each                |
| Need More?     | <a href="#">Request Quote</a> |

**i** Prices shown are exclusive of VAT/local taxes

Product Downloads

| General                          |               |
|----------------------------------|---------------|
| U-TLUIR-1-2 / 3-U840IEC          | Model Number: |
| Tube Lens                        | Type:         |
| Olympus                          | Manufacturer: |
| Physical & Mechanical Properties |               |

|                       |                             |
|-----------------------|-----------------------------|
| 350                   | Weight (g):                 |
| Optical Properties    |                             |
| 180.00                | Focal Length FL (mm):       |
| 1X                    | Magnification:              |
| 400 - 1000            | Wavelength Range (nm):      |
| 22.0                  | Field Number (mm):          |
| Regulatory Compliance |                             |
| View                  | Certificate of Conformance: |

## Product Details

- Empty Filter Cube and [Fluorescence Filter Pre-Mounted Version](#) Available
  - Wide Field and Single Port Tube Lenses
  - Low Auto Fluorescence Immersion Oil for High NA [Objectives](#)
  - C-Mount, T-Mount Adapters and Camera Coupler for System Integration
- Olympus Tube Lenses and Accessories are specifically designed for use with Olympus microscopes or objectives.


The Olympus Empty Filter Cubes can accept 25mm diameter [excitation and emission](#) filters and 25.2 x 35.6mm [dichroic filters](#) for customized fluorescence configurations. The [#86-832](#) Filter Cube is compatible with Olympus BX2 (BX61/51/41) and IX2 (IX81/71/51) series microscopes. The [#86-833](#) Filter Cube is compatible with Olympus BX3 (BX63/53/43) and IX3 (IX83/73/53) series microscopes. [Pre-Mounted Fluorescence Filter Cube Sets](#) utilizing the #86-833 Filter Cube are available for a variety of standard fluorophores.

Olympus Tube Lenses feature a 180mm focal length. The Olympus Wide Field Tube Lens ([#36-401](#)) offers a wide 26.5mm diameter image circle, and can be easily integrated into T-Mount systems with the Olympus Wide Field Tube Lens to T-Mount Holder ([#11-149](#))\*. The Olympus Single Port Tube Lens ([#86-835](#)) features a smaller 22mm diameter image circle, but is designed to easily integrate with microscope observation tubes for video observation. A C-Mount video microscope can be configured by using the #86-835 tube lens, [#89-606](#) Olympus C-Mount Adapter, and [#89-605](#) Olympus 1X Camera Coupler.

Olympus Low Auto Fluorescence Immersion Oil [#86-834](#) is designed for use with Oil Immersion Objectives to achieve a high numerical aperture. The immersion oil is specially formulated for reducing autofluorescence that is typically associated with oils, making it an ideal choice for Fluorescence Microscopy applications.

**\*Note:** The four socket head cap screws must be removed to mount the #36-401 Olympus Wide Field Tube Lens in the #11-149 Olympus Wide Field Tube Lens to T-Mount Holder. See Technical Information Tab.

## Technical Information



### Using Olympus Objectives with C-Mount Cameras

