

[See all 10 Products in Family](#)

Olympus SLMPLN 20X Objective

See More by [Olympus](#)



#86-824

Stock **#86-824** [CONTACT US](#)

⊖ 1 ⊕ €3.655⁰⁰

ADD TO CART

Volume Pricing	
Qty 1+	€3.655,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General

1-U2M265 **Model Number:**

Compatible Tube Lens Focal Length (mm):
Focal Length: 180mm

Microscope Objective **Type:**

Infinity Corrected **Style:**

Olympus
Manufacturer:

Physical & Mechanical Properties

Field of View (mm):
1.3

Length excluding Threads (mm):
20.00

Maximum Diameter (mm):
26

Weight (g):
56

Optical Properties

Compatible Cover Glass Thickness (mm):
N/A

Focal Length FL (mm):
9.00

Magnification:
20X

Numerical Aperture NA:
0.25

Transmission (%):
T > 90 from 400-700nm

Resolving Power (μm):
1.34

Depth of Field (μm):
4.40

Working Distance (mm):
25.0

Wavelength Range (nm):
400 - 700

Field Number (mm):
26.5

Parfocal Length (mm):
45

Immersion Liquid:
N/A

Threading & Mounting

Mounting Threads:
RMS / 20.32mm x 36 TPI

Regulatory Compliance


Certificate of Conformance:
[View](#)

Product Details


- Ideal for Brightfield Microscopy
- Ultra-Long Working Distances Reduce Risk of Specimen Damage
- Over 90% Transmission from 400 – 700nm
- Additional [Olympus Infinity Corrected Objectives](#) Available

Olympus Super Long WD MPlan Achromat Infinity Corrected Objectives are ideal for electronic assembly, semiconductor inspection, or other flat field inspection applications where resolution is critical. Olympus Super Long WD M Plan Achromat Infinity Corrected Objectives provide excellent performance, ultra-long working distances, and superior image clarity.


Technical Information



Typical Transmission



Using Olympus Objectives with C-Mount Cameras



Stock #ABCDEFG

Model	Part Number	20.0	0.5	0.2	16.1	25.0	26.0	4.9
SLMPLN 20X	86-824							
SLMPLN 50X	86-825							
SLMPLN 100X	86-826							

