

[See all 5 Products in Family](#)

Replacement Base for Peak Direct Measuring Microscopes



#52-822: Replacement Base for Peak Direct Measuring Microscopes

Stock **#52-822** **1 In Stock**

⊖ 1 ⊕ €32²⁵

ADD TO CART

Volume Pricing

| | |
|------------|-------------------------------|
| Qty 1-4 | €32,25 each |
| Qty 5+ | €30,70 each |
| Need More? | Request Quote |

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

Regulatory Compliance

Certificate of Conformance:

[View](#)

Product Details

- Multiple Magnification Options

- Illuminated Models Available with English or Metric Reticles
- Removable Clear Acrylic Base

Peak Direct Measuring Microscopes are high power, portable microscopes designed for precise depth measurements. Peak Direct Measuring Microscopes feature a 10X focusable eyepiece and a scale with an accuracy of 0.1mm for measuring the length, width, or depth of an object. The clear cutaway acrylic base allows for manipulation of the sample while viewing. Illuminated models feature a removable pen illuminator, and include two AA batteries. Peak Direct Measuring Microscopes are sold with a hard-shell vinyl storage case.

Note: Reticle should not be removed or replaced.

;