

In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union



SUBMIT

TECHSPEC® Shaft Collar (2 pk)



Shaft Collar (2 pk), #86-845

Stock #86-845 3-4 DAYS



1



€14,00

ADD TO CART

Qty 1+

€14,00

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Regulatory Compliance

RoHS:
Not Compliant

Certificate of Conformance:
[View](#)