

TECHSPEC[®] SpectrumOrange Filter Cube Set - Olympus



Olympus Filter Cube #86-833



Stock **#86-379-OLY** 1 In Stock

-

1

+

€1.350⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-9	€1.350,00 each
Qty 10+	€1.215,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General	
Fluorescence Filter Kit	Type:
SpectrumOrange	Compatible Fluorophore:
#67-032	Excitation Filter:

#84-101	Emission Filter:
#67-082	Dichroic Filter:
Olympus	Compatible Microscope Brand:

Optical Properties

Fused Silica (Corning 7980)	Substrate: ☐
Hard Coated	Coating:
532 - 554	Excitation Wavelength (nm):
576 - 596	Emission Wavelength (nm):
562.00	Dichroic Cut-On Wavelength (nm):
532 - 596	Wavelength Range (nm):

Regulatory Compliance

View	Certificate of Conformance:
------	-----------------------------

Product Details

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).