

TECHSPEC®

Texas Red™ Filter Cube Set - Olympus



Olympus Filter Cube #86-833



Stock **#67-008-OLY**

CONTACT US

-

1

+

€1.230^{,85}

ADD TO CART

Qty 1-9
€1.230,85
Qty 10+
€1.119,61
Volume Pricing
Request Quote
ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

SPECIFICATIONS

General

Fluorescence Filter Kit	
Type:	
Texas Red	
Compatible Fluorophore:	
Excitation Filter:	
#67-033	
Emission Filter:	
#67-035	
Dichroic Filter:	
#67-083	
Compatible Microscope Brand:	
Olympus	
Optical Properties	
Substrate: ☐	
Fused Silica (Coming 7980)	
Coating:	
Hard Coated	
Excitation Wavelength (nm):	
542 - 582	
Emission Wavelength (nm):	
604 - 644	
Dichroic Cut-On Wavelength (nm):	
593.00	
Wavelength Range (nm):	
542 - 644	
Regulatory Compliance	
RoHS:	
Not Compliant	
Certificate of Conformance:	
View	

PRODUCT DETAILS

- Ensure Proper Alignment of Filter Sets
 - Easily Adapt to Common Microscopes
 - Available in a Variety of Wavelengths
 - Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available
- TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.
- Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).
- Built-to-Order Product:** Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).