


[See all 1 Products in Family](#)

## UV Light Detection Label, 20/Pack



UV Light Detection Labels

Stock **#16-787** **5 In Stock**

1  €25<sup>00</sup>

**ADD TO CART**

### Volume Pricing

Qty 1-4	€25,00 each
Qty 5-9	€22,50 each
Need More?	<a href="#">Request Quote</a>

**!** Prices shown are exclusive of VAT/local taxes

### Product Downloads

#### General

Adhesive backing

**Note:**

#### Physical & Mechanical Properties

12.7 x 12.7

**Dimensions (mm):**

#### Optical Properties

## Regulatory Compliance

[View](#)

Certificate of Conformance:

[Compliant](#)

Reach 242:

## Product Details

- Self-Adhesive UV Light Detection Labels
- Permanent Record of UV Light Exposure
- Ideal for UV Disinfection, Curing, and Detection Applications

UV Light Detection Labels are self-adhesive labels that change color irreversibly upon exposure to UV light. The color density of the labels changes from white to blue with increased exposure time and distance, creating a permanent record of UV light exposure. UV Light Detection Labels are available in packs of 20 with 1/2 inch square dimensions and are ideal for germicidal disinfection, surface sterilization, UV curing, or general UV light detection applications.

## Technical Information

